

## ■ Specification

Performance	Trap Cooling Temperature	- 45°C
	Amount of Dehumidification	1 ℓ /one time
Function	Safety Function	Earth Leakage, Overcurrent Protection Breaker, Self-diagnosis Function, Refrigerator Protection (high-pressure pressure switch, overload relay), Vacuum Automatic Leak Valve, Service Outlet Fuse
Outer Dimensions/Weight		450 × 420 × 440mm / approx. 35kg
Power Input, Rated Power		6A 600VA · AC100V 50/60Hz

※ The value in performance is at 20°C as room temperature, rated supply voltage 50Hz.

※ Power input value does not include the capacity of service outlets.



콘텐츠산업 진흥법에  
의한 표시

1. 콘텐츠의 명칭 : 상품정보
2. 콘텐츠의 제작 및 표시 : 2017년 12월 22일
3. 콘텐츠의 제작자 : (주)라보텍
4. 콘텐츠의 이용조건
5. 콘텐츠의 보호기간 : (주)라보텍 사이트의 모든 콘텐츠는 [콘텐츠산업 진흥법]에 따라 제작 또는 갱신일로부터 5년간 보호됨

Copyright © 2014 by LABOTEC. All contents can not be copied without permission.